

Now everyone can connect to the conversation

Signia BCT IX



Universal
Connectivity



Enhanced User
Experience



Revolutionary
Sound Processing



signia



signia-pro.com



Signia BCT IX



Powered by our custom-designed BCT chip, the new Signia BCT IX is the perfect choice for you and your patients

Universal Connectivity	Enhanced User Experience	Revolutionary Sound Processing
Easily connects to 99% of devices including cell phones, tablets, laptops and TVs	Industry leading up to 36 hours of runtime including 5 hours of streaming and up to 54 hours without streaming	Track and enhance multiple talkers with our revolutionary RealTime Conversation Enhancement technology
Stream music, calls or entertainment straight to your device	The smallest Bluetooth Classic-compatible RIC hearing aid available today*	Ultra-fast, always-on e2e wireless 4.0 allows hearing aids to track and enhance multiple speakers in noisy environments
Stay connected on the go and answer/end HandsFree calls directly from your hearing aid	LED indicator light offers power status at a glance	OVP 2.0 ensures the wearers own voice is comfortable and provides a more pleasant and stable sound of both their own voice and surroundings
Integrated telecoil for better hearing in challenging environments supported by a loop system	Dual push buttons for volume and program control	Signia Assistant is the only live AI-driven companion that empowers the wearer to manage their own streaming and sound quality preferences
	Convenient visibility window to easily identify serial number and performance level	

Performance levels

7 IX

5 IX

3 IX

CROS IX

Choose from a range of colors



Key features

- Bluetooth Classic
- Longest li-ion runtime**
- IX Integrated Xperience
- RealTime Conversation Enhancement
- NFMI e2e wireless 4.0
- Own Voice Processing 2.0
- Signia Assistant - AI live DNN
- LED light
- Dual push buttons
- Telecoil
- IP68
- HandsFree calls
- Multiple charger options
- Recharge on the go

*In a comparison of RIC hearing aids based on the manufacturers' newest platforms as of 2024.

**Up to 54 hours without streaming. Up to 43 hours after 3 years of normal use.

Signia is a registered trademark of WSAUD A/S. SI-21683/0225

Be
Brilliant™

Unleash the Power of Conversation

Signia Integrated Xperience (IX)
with RealTime Conversation Enhancement



signia

Group Conversation in Noise is Difficult...

- Group conversations in noise are the biggest unmet need for people with hearing loss
- Can lead to frustration, isolation and ultimately withdrawal
- They miss out on the moment, we miss out on their contribution to the conversation



Signia's Integrated Xperience provides a revolutionary three-stage process that puts conversation first.

ANALYZE



Evaluates 192,000 data points every second to build a picture of the listening situation and to understand the environment and conversation needs

AUGMENT




Boosts the voices of all the people in the conversation, helping the wearer to hear them clearly while staying immersed in their surroundings

ADAPT



Updates 1,000 times per second, allowing it to seamlessly adjust in real time as the conversation moves and shifts, ensuring the wearer is always part of the action



Signia Integrated Xperience boasts an impressive 3.2 dB SNR improvement over the nearest competitor in a noisy group conversation scenario.*

*Visit signia-pro.com/en-us/claims-data for more details.

Benefits of RealTime Conversation Enhancement*:

A comfortable,
immersive
conversational
experience without
distractions

Access to multiple
speakers via the
world's first multi-
stream architecture¹
with speech that
pops above the
surrounding sound

Keep up and
contribute to the
conversation

*RealTime Conversation Enhancement is not available in all styles and performance levels.


1. Number of focus streams varies per technology level.

2. Visit signia-pro.com/en-us/claims-data for more details



95%

of wearers have better
speech reception with
RealTime Conversation
Enhancement than without.²



Signia BCT IX

The perfect choice for your patients

What makes this the perfect choice for you and your patients? Simple and effective hearing care thanks to universal connectivity, an enhanced user experience, and revolutionary sound processing.

Universal connectivity

- Easily connects to 99% of devices including cell phones, tablets, laptops & TVs
- High quality audio-experience streaming music, calls and entertainment
- Stay connected on the go and answer/end HandsFree calls directly from the hearing aid
- Integrated telecoil for better hearing in challenging environments supported by a loop system



Revolutionary Speech Processing

- Our RealTime Conversation Enhancement delivers 22% better speech understanding in noisy group conversations than leading competitors**
- The only Bluetooth classic compatible hearing aid with ultra-fast, always-on NFMI e2e wireless 4.0 allows hearing aids to track and enhance multiple speakers in noisy environments.

Enhanced user experience

- Industry leading 36 hours of runtime — including 5 hours of streaming; to 54 hours without***
- Convenient visibility window to easily identify the serial number and performance level

Available Chargers:



Available Colors:



Available Performance Levels: 7IX, 5IX, 3IX, compatible with CROS Pure BCT IX

*In a comparison of RIC hearing aids based on the manufacturers' newest platforms

** Jensen NS, Mohnlein-Gilbert K, Wilson C, Berwick N, Kamkar Parsi H, Best S, Samra B & Taylor B. 2024. Signia IX with pioneering multi-stream technology delivers 22% better speech understanding in noisy group conversation than a competitor with an AI co-processor-driven platform. Signia White Paper.

*** Up to 54 hours without streaming

Silk Charge&Go IX

The only thing that changes is your hearing

Silk Charge&Go IX – an innovative completely-in-canal (CIC) hearing aid that pairs exceptional sound quality with unmatched discretion.



Power through the day

The convenience of a rechargeable battery for all-day listening. Plus, the accompanying hassle-free case charges wirelessly.

Binaural OneMic Directionality

The only CIC device that offers Binaural OneMic Directionality to focus on a speaker in noisy listening environments.

Battery run time up to 28 hours

Conveniently charge hearing aids up to 4 times on the go.

Unmatched discretion

The soft, textured surface avoids reflections for added discretion.

Immediate fit, immediate impact

The quick, satisfying experience of a ready-to-wear CIC. The Silk Charge&Go IX's small, snug design is discreet and it can be worn immediately.

Control with the app

The user-friendly Signia app and miniPocket remote control allows patients to fine-tune their hearing aids to their unique needs for personalized control anytime, anywhere.

Available Faceplate Colors:



Black



Mocha

Available Performance Levels: 7IX, 5IX, 3IX, 2IX, 1IX, compatible with CROS Silk Charge&Go IX

Silk Charge&Go IX has an IP68 rating ensuring the hearing aids will still function if immersed under 1 m of water for 1 h

Insio IX IIC/CIC

Filter out the noise – keep the conversation

Our classic Insio CICs and IICs just got better! Now with a new design and built on our latest platform. Insio IX IIC/CIC is the only custom one-mic hearing aid offering binaural directionality for focused conversation, even in noisy environments.

Advanced clarity

Newly upgraded and meticulously crafted to comfortably enhance speech clarity while in background noise.

Modern tech

Our Binaural OneMic Directionality 2.0 helps focus on speech in noise.

Ultimate discretion

Our smallest custom hearing aids, available in both Completely-in-Canal (CIC) and Invisible-in-Canal (IIC) models.

Enhanced comfort

The new Insio IX IIC/CIC models come with an improved design for a snug, comfortable fit.

Available Colors:



Transparent Red
(shell)/Mocha
(faceplate)



Transparent Red
(shell)/Deep
Brown (faceplate)



Transparent Red
(shell)/Black
(faceplate)



Transparent Blue
(shell)/Mocha
(faceplate)



Transparent Blue
(shell)/Deep
Brown (faceplate)



Transparent Blue
(shell)/Black
(faceplate)

Available Performance Levels: 7IX, 5IX, 3IX, 2IX, 1IX, compatible with CROS Silk Charge&Go IX

Styletto IX



A conversation starter

Styletto IX keeps the conversation flowing without compromising on style. Available in various colors, this sleek hearing aid is the refined choice for better hearing.

Signia Integrated Xperience

Designed to adapt to the wearer's conversational environment in real time, so they can take part with confidence.

Power through the day

The convenience of a rechargeable battery that supports hearing all day with battery runtime up to 20 hours. Plus, the portable charging case holds up to 5 full charges, enabling on-the-go charging.

Expanded connectivity

- Bluetooth LE Audio
- ASHA (Audio Streaming for Hearing Aids)
- Made for iPhone®
- Bluetooth® connectivity via streamer

Nine distinctive colors and an ultra-sleek design.

Personalized control

Easy-to-use mobile app lets wearers fine-tune their hearing aids whenever and wherever needed.

Available Chargers:



Styletto Portable Charger



Dry&Clean Charger

Available Colors:



Available Performance Levels: 7IX, 5IX, 3IX, 2IX, 1IX, compatible with CROS Styletto IX

Styletto IX has an IP68 rating ensuring the hearing aids will still function if immersed under 1 m of water for 1h (excluding receiver).

Active Pro IX / Active IX

Reconnect with the conversation

Welcome to a new era of hearing with Active Pro IX — a game-changing prescription hearing aid featured in a modern earbud design. Includes RealTime Conversation Enhancement* and portable rechargeability. Fits instantly, enables hands-free calls and streaming, and resists sweat and dust to support an active lifestyle.

All-day power

Wearers receive up to 34 hours of wear time and for added convenience the Qi wireless portable charging case can provide up to 3 full charges on the go.

Bluetooth connectivity

Can be paired to compatible Apple® iOS, Bluetooth LE Audio, and ASHA-supported Android™ devices for direct streaming of audio such as phone calls and music.




Water & dust resistant

Active Pro IX and Active IX are IP68 certified so wearer can enjoy an active life without worrying about sweat, water, or dust.

HandsFree phone calls with CallControl

Patients can answer or end phone calls with a simple tap of the hearing aid enabling HandsFree calling for compatible MFi and Bluetooth LE Audio enabled devices.

Available Colors:



White

Rose Gold

Black

Brown

Mocha

Beige

Available Performance Levels: Active Pro IX and Active IX, compatible with CROS Active IX

* RealTime Conversation Enhancement only available in Active Pro IX.

Active Pro IX, Active IX, and CROS Active IX have an IP68 rating ensuring the hearing aids will still function if immersed under 1 m of water for 1h.



Future-proof connectivity & performance



With Signia, you give your patients the chance to stay connected with the future of hearing care.

- Bluetooth Classic
- Bluetooth LE Audio available for Bluetooth-enabled Integrated Xperience hearing aids
- ASHA (Audio Streaming for Hearing Aids)
- Made for iPhone®
- Bluetooth connectivity via streamer

Signia also provides firmware updates to enhance performance.

Own Voice Processing 2.0*



Own Voice Processing (OVP) promotes hearing wellness

Contribute comfortably — OVP 2.0 ensures the wearers own voice is comfortable and provides a more pleasant and stable sound of both their own voice and surroundings.

Boost Communication — 78% of participants reported increased communication when using OVP⁵

Signia Assistant*



The 24/7 hearing companion

Signia Assistant is the only live AI-driven companion that empowers the wearer to manage their own streaming and sound quality preferences.



* Not available in all styles and performance levels.

5 Heydal EH. A new own voice processing system for optimizing communication. Hearing Review.



Learn more today.



"Made for iPhone" mean that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone may affect wireless performance. iPad is a trademark of Apple Inc., registered in the U.S. and other countries.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc.

signia-pro.com

Signia is a registered trademark of WSAUD A/S. SI-21711

Be
Brilliant™

Now everyone can connect to the conversation



— with the groundbreaking
Signia Pure Charge&Go
BCT IX



signia

Pure Charge&Go BCT IX



Join the conversation



Dive into the world of enhanced hearing with the Pure Charge&Go BCT IX.

Scan the QR code now to uncover all its features — or speak directly with your Hearing Care Professional.

Don't wait! Get brilliant hearing today with the Pure Charge&Go BCT IX.



Discover the perfect blend of comfort, stylish design, and innovative technology with Signia.



Connect your way

It's the perfect choice for unmatched hearing clarity and seamless connectivity.





Why choose the BCT IX?


- **Stay connected** on the go and answer/end **HandsFree** calls directly from your hearing aid.
- With an unprecedented **36 hours of battery life**, including **5 hours of streaming**, you have the freedom to go about your day.
- The **smallest Bluetooth Classic-compatible RIC*** hearing aid features a lightweight, sleek design that even the most discerning wearers will appreciate.
- The **new LED light** instantly shows your hearing aids' power status.
- **Dual push buttons** allow quick volume and program adjustments on the go.





The new Signia Pure Charge&Go BCT IX keeps you confidently connected to the conversation

Powered by the Signia IX platform, the groundbreaking RealTime Conversation Enhancement delivers **22% better speech understanding** in noisy group conversations than leading competitors.¹ So you can contribute to the conversation, no matter how busy it gets.



About Signia

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath. With Signia, you can Be Brilliant.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners. Other trademarks and trade names are those of their respective owners.

¹ Jensen NS, Mohnlein-Gilbert K, Wilson C, Berwick N, Kamkar Parsi H, Best S, Samra B & Taylor B. 2024. Signia IX with pioneering multi-stream technology delivers 22% better speech understanding in noisy group conversation than a competitor with an AI co-processor-driven platform. Signia White Paper.

* In a comparison of RIC hearing aids based on the manufacturers' newest platforms.

Signia is a registered trademark of WSAUD A/S. SI-21710/0225

Find us on:



signia.net



facebook.com/signiahearing



youtube.com/signiahearing



instagram.com/signiahearing